

HCP-7000
CHART PROJECTOR

HUVITZ
CHART PROJECTOR



HCP-7000 CHART PROJECTOR

- Semi-permanent lifetime LED light source
- Brighter and clearer than conventional bulb light
- Attractive Streamlined Shape and Blue Indicator
- RC-R Type (Russian Chart) added



SPECIFICATIONS

Chart	41 charts 34 masks, Red/Green & Polarization Filters
Projection Distance	2.5~8m
Chart Rotation Speed	Average 0.15 sec
Projection Magnification	30x at 5m
Power Saving	Automatic Switch off (10 min)
Program	2 programs with a maximum of 30 charts each
Tilt Angle	15°
Power Supply	100-120V 50Hz : 0.6A, 200-240V 60Hz : 0.3A
Lamp	LED 4W
Dimension	270(W) x 182(D) x 230(H)mm
Weight	3.44kg

Designs and details can be changed without prior notice for the purposes of improvement.



HUVITZ Co., Ltd. 38, Burim-ro 170beon-gil, Dongan-gu,
Anyang-si, Gyeonggi-do, 14055, Republic of Korea
Tel:+82-31-442-8868 Fax:+82-31-477-8617 <http://www.huvitz.com>

E2ARAA-16-00001-2-170426