

# HUVITZ APPLANATION TONOMETER HT-5000





# HUVITZ APPLANATION TONOMETER HT-5000

- The world standard Goldmann applanation tonometer
- The digital display makes it easy to read even in a dark room.
- The wide grip comfortable knob for easier adjustment of measuring prism
- Optimized algorithm realizes the fast process of IOP calculation.
- Instantaneous wake-up time from sleep mode!
- The power supply by the most commonly used AAA batteries guarantees easy installation and maintenance.
- The slanted display matches to the angle of a user's view for more comfortable observation.

## SPECIFICATIONS

Measuring force	generated by leverage weight
Installation	Mountable on peg on microscope.
Measuring range	3~75mm Hg
Reverse span	$\pm < 0.49$ mN
Power Consumption	DC3V
Net weight	0.50kg (without accessories, with batteries)

Designs and details can be changed without prior notice for the purposes of improvement.

**Huvitz**

HUVITZ Co., Ltd. 38, Burim-ro 170beon-gil, Dongan-gu,  
Anyang-si, Gyeonggi-do, 14055, Republic of Korea  
Tel: +82-31-442-8868 Fax: +82-31-477-8617 <http://www.huvitz.com>

M1ARAA-16-00001-2-170426